

# Wide-neck bottle and Multipurpose containers



## Wide-neck bottle, square bottom

High quality polyethylene (PE-HD).

- color: natural, translucent
- for easy identification the closure is made of blue polypropylene
- with originality seal ring on the closure
- 250 ml and 500 ml bottle with UN approval 3H1/x0.8/S
- supplied: bottle and separate closure, non-sterile

Volume ml	Fig. No.	Neck Ø mm	Height mm	W x D mm	Packaging	Quantity per Pack	Order No.
50	1	28	88	36 x 36	cardboard box 1 x 100	100	59 24 050
250	2	45	129	59 x 59	cardboard box 1 x 100	100	59 24 250
500	3	54	162	75 x 75	cardboard box 1 x 50	50	59 24 500

**New**  
**New**  
**New**



## Multipurpose containers

Polypropylene (PP).

The best round containers for sampling or storing of liquids, paste or powder medicine. The containers are available in three sizes.

- for easy identification and the option to mark them with a pen the containers are white
- suitable for freezer storage
- open autoclavable at 121 °C
- supplied: container with screw-cap, non-sterile

Volume ml	Fig. No.	Ø mm	Height mm	Packaging	Quantity per Pack	Order No.
50	1	35	74	cardboard box 1 x 200	200	59 23 050
250	2	65	90	cardboard box 1 x 100	100	59 23 250
500	3	90	86	cardboard box 1 x 72	72	59 23 500

**New**  
**New**  
**New**