

ratiolab® Cleanibag System and Autoclavable bags BIOHAZARD



ratiolab® Cleanibag System

All-purpose PE bags with self-adhesive strip for clinical and laboratory use.

- self-adhesive strip, sticks on smooth and clean surfaces (e.g. walls, work benches, beds, bedside tables)
- residue-free on removable

On request: other sizes, including for intensive care, bags with custom-made imprint

W x H mm	Volume approx. l	Packaging	Quantity per Pack	Order No.
200 x 300	1.5	cardboard box 1 x 1000	1000	70 00 250
300 x 300	6.5	cardboard box 1 x 1000	1000	70 00 260
300 x 450	8	cardboard box 1 x 1000	1000	70 00 270

Indicator patch



ratiolab® Autoclavable bags BIOHAZARD red – Warning Instruction and sterilization indicator patch

Made of robust polyethylene (PE).

Distinctive red autoclavable bags for easy identification, for the disposal of contaminated labware.

- rectangle white indicator field turns into dark when autoclaved
- with BIOHAZARD symbol and safety imprint in English, German, French and Spanish languages
- autoclavable at 121 °C, for autoclaving, the bag must not be completely closed
- practical bench-top stand, plastic-coated, with 3 rubber feet



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
210 x 290	1.5	20 µm	cardboard box 10 x 100	1000	New 70 02 110
Set consisting of bench-top stand and 100 red waste disposal bags 210 x 290 mm				1	New 70 01 410

ratiolab® Autoclavable bags BIOHAZARD – Warning Instruction and sterilization indicator patch

Autoclavable disposable waste bag of very strong PP-film with a broad seam at the bottom, for the sterilization and disposal of hazardous and contaminated labware.

- rectangle white indicator field turns into dark when autoclaved
- with BIOHAZARD sticker and safety imprint in English, German, French and Spanish language
- autoclavable at 134 °C, for autoclaving, the bag must not be completely closed



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
300 x 500	8	50 µm	cardboard box 1 x 500	500	70 02 100
400 x 780	38	50 µm	cardboard box 1 x 500	500	70 02 150
600 x 800	60	50 µm	cardboard box 1 x 500	500	70 02 200
700 x 1100	110	50 µm	cardboard box 1 x 350	350	70 02 300

ratiolab® Waste Disposal Bag

ratiolab® Waste Disposal Bag, Standard

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 134 °C, for autoclaving, the bag must not be completely closed
- only for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet



W x H mm	Volume approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
200 x 300	2.5	50 µm	dispenser 1 x 100	100	70 01 000
300 x 500	8	50 µm	cardboard box 1 x 500	500	70 01 100
400 x 780	38	50 µm	cardboard box 1 x 500	500	70 01 150
600 x 800	60	50 µm	cardboard box 1 x 500	500	70 01 200
600 x 800	60	50 µm	cardboard box 1 x 100	100	70 01 210
600 x 800	60	100 µm	cardboard box 1 x 125	125	70 01 212
700 x 1100	110	50 µm	cardboard box 1 x 350	350	70 01 300
700 x 1100	110	50 µm	cardboard box 1 x 75	75	70 01 310
700 x 1100	110	100 µm	cardboard box 1 x 75	75	70 01 311
bench-top stand for waste disposal bags 200 x 300 mm				1	70 01 400
Set consisting of bench-top stand and 100 waste disposal bags (standard) 200 x 300 mm				1	70 01 401

ratiolab® Waste Disposal Bag, High Transparency

Bags of very strong PP-film with a broad seam at the bottom, for the disposal of contaminated labware.

- autoclavable at 121 °C, for autoclaving, the bag must not be completely closed
- only for bag size 200 x 300 mm: practical bench-top stand, plastic-coated, with 3 rubber feet



W x H mm	Volumen approx. l	Film Thickness	Packaging	Quantity per Pack	Order No.
200 x 300	2.5	50 µm	dispenser 1 x 100	100	70 01 005
300 x 500	8	50 µm	cardboard box 1 x 500	500	70 01 105
400 x 780	38	50 µm	cardboard box 1 x 500	500	70 01 155
600 x 800	60	50 µm	cardboard box 1 x 500	500	70 01 205
700 x 1100	110	50 µm	cardboard box 1 x 350	350	70 01 305
bench-top stand for waste disposal bags 200 x 300 mm				1	70 01 400
Set consisting of bench-top stand and 100 waste disposal bags (high transparency) 200 x 300 mm				1	70 01 405