

Multi Purpose Beakers

Multi Purpose beakers with screw-caps



Polypropylene (PP) or polystyrene (PS).

For sampling, transport and storage of solid, pasty and liquid samples.

- with red screw-cap
- supplied in non-sterile or sterile, individually packed



Volume ml	Fig. No.	Material	Graduation and Labeling area	Sterilized	Ø* x Height mm	Packaging	Quantity per Pack	Order No.
30	1	PP	–	–	49 x 43	cardboard box 1 x 1000	1000	59 21 173
30	1	PP	–	■	49 x 43	bags 1000 x 1	1000	59 21 174
60	2	PP	■	–	67 x 43	cardboard box 1 x 600	600	59 21 175
60	2	PP	■	■	67 x 43	bags 500 x 1	500	59 21 176
120	3	PP	■	–	73 x 62	cardboard box 1 x 500	500	59 21 182
120	3	PP	■	■	73 x 62	bags 250 x 1	250	59 21 183
200	4	PP	–	–	82 x 64	cardboard box 1 x 200	200	59 21 185
200	4	PP	–	■	82 x 64	bags 150 x 1	150	59 21 184
200	5	PS	–	■	82 x 64	bags 150 x 1	150	59 21 188

*with screw-cap

ratiolab® Multi Purpose Beaker System and Multi Purpose bottles



ratiolab® Multi Purpose Beakers

Glass-clear polystyrene (PS).

- volume 100 ml, with graduation
- self-adhesive film cover, leakproof and easily labelled
- comfortable sampling, the film cover can be pierced with pipet tips
- practical packaging, 50 multi purpose beakers are stacked in PE-bags for use with the beaker dispenser

Product	Packaging	Quantity per Pack	Order No.
ratiolab® multi-purpose beakers with cover	bags 20 x 50	1000	59 21 170
ratiolab® multi-purpose beakers without cover	bags 20 x 50	1000	59 21 171
self-adhesive film covers for beakers	bag 1 x 500	500	59 21 160

ratiolab® Beaker Tray and Beaker Dispenser for Multi Purpose Beakers

Beaker Tray: Robust plastics with easy-to-clean surface.

- for 10 ratiolab® multi purpose beakers

Beaker Dispenser: Robust plastics with easy-to-clean surface.

- for 50 ratiolab® multi purpose beakers
- for individual beaker dispensing
- glass-clear dispenser door for refills makes stock visible
- separate side box for stocking of self-adhesive film covers

Product	Fig. No.	Quantity per Pack	Order No.
ratiolab® beaker tray	1	1	59 21 200
ratiolab® beaker dispenser	2	1	59 21 190

Multi Purpose bottles

Polyethylene, blue translucent (PE).

- inclusive lid with nonwoven fabric liner
- with graduations, light protection
- wide neck diameter 80 mm
- Volume: 2.5 l

Product	Packaging	Quantity per Pack	Order No.
Multi purpose bottle 2.5 l with self-adhesive labels	cardboard box 1 x 60	60	59 22 250
lids for Multi purpose bottle 2 and 2.5 l	bag 1 x 100	100	59 22 204
self-adhesive labels for Multi purpose bottle 2.5 l	bag 1 x 100	100	59 22 253