

ratiolab® Multistand Rack and Rack for Microcentrifuge Tubes




ratiolab® Multistand Rack

Polyamide (PA).

With the ratiolab® Multistand you have the opportunity to use it with several microcentrifuge tubes, cryo tubes and test tubes with different diameters through the patented holder to fix safely.

- for tubes Ø 10 to Ø 17 mm
- dishwasher safe
- max. working temperature: 90 °C



The special designed inner clamps (patent) fixes the test tubes centric and vertical at any time.



Color	Positions	W x D x H mm	Packaging	Order No.
white	5 x 5	105 x 105 x 42	1	37 12 101



Rack for Microcentrifuge Tubes

Polypropylene (PP).

Ideal for the secure handling of reaction vessels. Space-saving on the laboratory bench due to its narrow form. Practical side grips.

- 5 x 16 wells for 1.5 and 2.0 ml reaction vessels
- autoclavable at 121 °C
- temperature-resistant from – 80 °C to +121 °C
- 225 mm x 65 mm x 28 mm (LxWxH)

Color	For micro tubes	Positions	Packaging	Order No.
green	1.5 – 2.0 ml	5 x 16	5	37 12 006
pink	1.5 – 2.0 ml	5 x 16	5	37 12 008
blue	1.5 – 2.0 ml	5 x 16	5	37 12 005
orange	1.5 – 2.0 ml	5 x 16	5	37 12 007
yellow	1.5 – 2.0 ml	5 x 16	5	37 12 009